



Dual

Disinfectant & Deodorant

Dual is for use on pre-cleaned, hard, non-porous, non-food contact surfaces in/on industrial, institutional, commercial, medical, and residential facilities and vehicles. Hospital use and broad-spectrum disinfectant and antimicrobial agent. Two-way deodorizing action. Pleasantly scented. Attacks malodors at their source. Sanitizes fabric and leather. Controls mold, mildew, and fungi. EPA Registered.

DIRECTIONS FOR USE

GENERAL: Shake container well before each use and hold it upright while spraying. For use on pre-cleaned, hard non-porous surfaces, not for use on humans or animals. Not recommended for use on waxed/polished surfaces, water-based latex paint, shellac, cellulose (such as rayon) r other alcohol-sensitive materials.

FOR DISINFECTING, MILDEW CONTROL, SANITIZING, & DEODORIZING OF CLEAN NON-FOOD CONTACT SURFACES: Hold container upright about 6 - 8" (1/5m) from surface and spray until target, area is thoroughly wet. Complete disinfection and tuberculocidal activity require a contact time of 10 minutes at approximately 68°F (20°C). Sanitization requires a 5-minute contact time at about 72°F (22°C). For mildew control, repeat application when necessary.

See label for further instructions.

TECHNICAL DATA

Color.....	Clear
Scent.....	Pleasant
Foaming.....	N/A
Solubility in Water	Complete
Density (lb/gal)	N/A
Flash Point	123°F
pH	N/A





FEATURES & BENEFITS

Dual

Disinfectant & Deodorant

FEATURES

Fresh lemon scent

Many uses

Tuberculocidal claim

Full 20 ounce can

CUSTOMER BENEFITS

No hospital smell

Disinfects and deodorizes in one step

Meets bloodborne pathogens rule

Our can last longer than competitions smaller cans

