

Dual

Disinfectant & Deodorant

Dual is for use on pre-cleaned, hard, non-porous, non-food contact surfaces in/on industrial, institutional, commercial, medical, and residential facilities and vehicles. Hospital use and broad-spectrum disinfectant and antimicrobial agent. Two-way deodorizing action. Pleasantly scented. Attacks malodors at their source. Sanitizes fabric and leather. Controls mold, mildew, and fungi.

Directions for Use

MOPPING GENERAL: Shake container well before each use and hold it upright while spraying. For use on pre-cleaned, hard nonporous surfaces, not for use on humans or animals. Not recommended for use on waxed/polished surfaces, water-based latex paint, shellac, cellulose (such as rayon) or other alcoholsensitive materials.

FOR DISINFECTING, MILDEW CONTROL, SANITIZING, & DEODORIZING OF CLEAN NON-FOOD CONTACT SURFACES: Hold container upright about 6 - 8" from surface and spray until target, area is thoroughly wet. Complete disinfection and tuberculocidal activity require a contact time of 10 minutes at approximately 68°F (20°C). Sanitization requires a 5-minute contact time at about 72°F (22°C). For mildew control, repeat application when necessary.

See label for further instructions.

Safety Reminder

Consult product label and Safety Data Sheet (SDS) before use. SDS's may be downloaded from www.gatorchemical.com.

> Gator Chemical 2202 Industrial Blvd Sarasota FL, 34234 800-234-3672 www.gatorchemical.com

Technical Specifications

	Color	Clear
	Scent	Pleasant
t	Foaming	N/A
-	Solubility in Water	Complete
	Density (lb/gal)	N/A
	Flash Point	123°F
	рН	N/A

Packaging

12 x 15 Ounces

