

Dual

Disinfectant & Deodorant

Dual is for use on pre-cleaned, hard, non-porous, non-food contact surfaces in/on industrial, institutional, commercial, medical, and residential facilities and vehicles. Hospital use and broad-spectrum disinfectant and antimicrobial agent. Two-way deodorizing action. Pleasantly scented. Attacks malodors at their source. Sanitizes fabric and leather. Controls mold, mildew, and fungi.

Directions for Use

MOPPING GENERAL: Shake container well before each use and hold it upright while spraying. For use on pre-cleaned, hard non-porous surfaces, not for use on humans or animals. Not recommended for use on waxed/polished surfaces, water-based latex paint, shellac, cellulose (such as rayon) or other alcohol-sensitive materials.

FOR DISINFECTING, MILDEW CONTROL, SANITIZING, & DEODORIZING OF CLEAN NON-FOOD CONTACT SURFACES: Hold container upright about 6 - 8" from surface and spray until target, area is thoroughly wet. Complete disinfection and tuberculocidal activity require a contact time of 10 minutes at approximately 68°F (20°C). Sanitization requires a 5-minute contact time at about 72°F (22°C). For mildew control, repeat application when necessary.

See label for further instructions.

Safety Reminder

Consult product label and Safety Data Sheet (SDS) before use. SDS's may be downloaded from www.gatorchemical.com.

Technical Specifications

Color	Clear
Scent	Pleasant
Foaming	N/A
Solubility in Water	Complete
Density (lb/gal)	N/A
Flash Point	123°F
pH	N/A

Packaging

12 x 15 Ounces



Gator Chemical
2202 Industrial Blvd
Sarasota FL, 34234
800-234-3672
www.gatorchemical.com

